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PRESENTACIÓN  
PRESENTATION  
PRESENTAZIONE

SOLUTIONS

# KITBRETT

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developed by

VIVELYYS  
WINE ECO-ENGINEERING

## SOMMAIRE

### SOMERO - CONTENTS - SOMMARIO

Présentation du KITBRETT .....	3-6
Presentación del KITBRETT .....	7-10
Presentation of the KITBRETT .....	11-14
Présentazione de KITBRETT .....	15-18

SOLUTIONS

# KITBRETT

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SOLUTIONS

# KITBRETT

Detection and prevention of wine contamination  
by *Brettanomyces*

## PUT YOUR WINE UNDER HIGH SURVEILLANCE

Since a few years, the consumer informed or not research **fruity characters in the wines which it tastes and becomes very attentive with their aromatic quality.** One regularly observes the rejection of the characters leather, inks, sweat of horse pronounced, which are the result of the growth in the wine of the *Brettanomyces* yeast.

In all the cases, **the activity of *Brettanomyces* leads to the loss of the fruity characters, a total draining of the wine in mouth and in its term, with the appearance of unpleasant odor.** Thus this yeast can in little time to destroy all the efforts of breeding of the producer and became a pet peeve for the oenologist.

The period between the end of alcoholic fermentation and the FML is critical for the appearance of contaminant because the wine is not protected by sulphur. *Brettanomyces* can benefit from this period to develop if the conditions are favorable.

After sulphiting, detrimental growths can also be observed if the quantity of active SO<sub>2</sub> is insufficient.

These yeasts develop then actively and the wine requires a treatment which will be function of the level of population and the stage of vinification.

**Only a systematic and regular analytical follow-up of the populations makes it possible to take the good decisions at the good moment, before the appearance of problems.**

FRANÇAIS

ESPAÑOL


ENGLISH

ITALIANO

## Description of **KITBRETT Complet**

«KITBRETT Complet» for *Brettanomyces* contamination's detection contains all the required equipment for 40 reliable analyses.

For those who are already equipped, KITBRETT Milieu contains 40 culture medium alone.

## Protocol of analysis:

The new formula of Vivelys medium culture allows **early and fast tracking of *Brettanomyces***, since the end of the fermentation with a reliable recognition thanks to the characteristic odor animal in case of positive analyse. This medium is also usable for the analysis of the ciders and beers.

One estimates necessary **5 to 8 analyses per batch and a year** for an effective protection of the wines according to the level of risk of the winery.

For a regular follow-up, it's important to carry out the first analysis of a batch as of the end of alcoholic fermentation for a maximum effectiveness of the follow-up. The test sample selection must be made in a sterile way to guarantee the reliability of the analyses.

1. Mark the name of the sample as well as the date of sampling on the box.
2. After homogenization, take 0.1 to 0.5 mL of the sample to be analyzed (wine, dregs...) using the sterile pipette (single use, 13 drops = 0.5 mL).
3. Open the box under a flame and spread out the wine over agar.
4. Place the box at the drying oven during 4 to 5 days at 28°C.

The culture medium can also be used for the quality controls of the wines bottled or the clear wines at the end of the breeding before setting. For that :

1. Filter vacuum in sterile environment, 10 to 100 mL of the wine to be analyzed, on nitrocellulose membrane 0.45 µm.
2. Deposit delicately the membrane on the surface of agar face impregnated upwards.
3. Incubate the box with 28°C during 4-5 days.

## Counting and interpretation of the results:

According to the stage of vinification, it is possible that certain stocks of another micro-organism different of *Brettanomyces* are able to develop on the milieu. It is currently the case of all the culture media *Brettanomyces* present on the market.



Vivelys KITBRETT new culture medium has a system of recognition in the event of presence of *Brettanomyces*: a precursor aromatic for an olfactive recognition (spiced odor or animal characteristic).

*Brettanomyces* colonies are characterized by their size, shape, edges, colour and the morphology of the yeast: **small white colonies, smooth and creamy** (photo 1), appearing late (between 4 or 5 days).

In the event of doubt, microscopic examination of the colonies (magnification 40x10) establishes recognition (photo 2).

The result expressed in UFC/mL (Unit Forming Colony) is given by the ratio: number of colony/volume of spread out sample.

For 0.5 mL of spread out wine, if 5 colonies are observed, one estimates that the cellular concentration of the starting sample is of 10 UFC/mL.

Photos 1 and 2: result examples after a culture during 4 days, 1) *Brettanomyces* colonies, 2) *Brettanomyces* cells on microscop.

When the number of colony exceeds the 100, countings become long and tiresome. It is then preferable to make an estimate. One can use the squaring printed on transparent paper delivered with KITBRETT.

This squaring is made of 20 large squares of which some are subdivided in 9 small squares.

Place the squaring on the bottom of the box. To count the number of colonies on 5 larges squares chosen randomly and make the average. For 0.5 ml of spread out wine, the result in UFC/mL will be Nb of colonies average x 2 x 20 = Nb of colonies average X 40.

**Caution:** Vivelys has produced this microbiological analysis kit for *Brettanomyces* to provide a simple and fast means of performing a check, but does not guarantee absolute detection of the organisms concerned. This product is used at the user's own risk. Vivelys cannot be held responsible, directly or indirectly, for the detection or non-detection of the organisms excluded by this test. Users must rule out the possibility of taking any action or pressing any charges against Vivelys for errors or damage caused to any product following detection or non-detection of the presence of suspect organisms.

Careful: the medium contains antibiotics, to avoid the contact with the skin, not to swallow.  
Conservation of KITBRETT Milieu: between 4 and 24°C in vacuum packed.  
To use before end: see on the side of the Petri dish.

**To introduce *Brettanomyces* populations monitoring in your winery, Vivelys also offers you training and annual consulting services: optimization of microbiological analyses, SO<sub>2</sub> management, optimization of winery hygiene, full audit especially.**

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